

Crutches

Model Number: 10401-8

\$34 Each

Walking Crutches with Underarm Pad and Handgrip



- Comfortable, durable underarm pad and hand grip
- Double extruded center tube provides additional strength to weight-bearing areas
- Easy wing nut adjustments for hand grips
- Lighter than wood, the standard aluminum crutch provides stability and durability
- Non-skid, jumbo size vinyl contoured tips provide excellent traction
- Push-pin adjustment makes it easy to adjust crutch height in 1" increments

❖ Cane Length: 45.75"

For Additional Information Please Click [Here](#)

Crutches

Model Number: 10403HD

\$125

Lightweight Walking Forearm Crutches



- Extra large tips for added stability
- Leg and forearm sections adjust independently for optimal sizing
- Vinyl coated, ergonomically contoured arm cuffs molded for comfort and stability
- Vinyl hand grips are comfortable and durable

For Additional Information Please Click [Here](#)

Weight (of the item)	Weight Limit
5.7 lbs	500 lbs

Crutches

Model Number: 10410

\$51

Euro Style Light Weight Forearm Walking Crutch



- Available 2 in attractive colors
- One piece molded plastic cuff and hand grip assembly provides safety and comfort

For Additional Information Please Click [Here](#)

Weight (of the item)	Weight Limit
1.5 lbs	300 lbs