

**Power Chairs**

**Model #: GP-162**

**\$1,425**

✓ Free Delivery \*



- Unique, under seat storage basket
- Comfortable, full-size stadium style seat
- Black, non-scuffing tires with shiny hubcaps on steel rims
- Smooth 26" turning radius great for indoors
- Comes standard with 2 22AH batteries
- Operating range up to 15.5 miles
- Top speed 3.5 mph

For Additional Information Please Click [Here](#)

**LiteRider Envy**

Weight Capacity	Seat Width/Depth	Ground Clearance	Driving Range	Max Speed
300 lbs	17"/16"	2.5"	15.5 miles	3.5 mph

\*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00  
(mileage charges may apply-see above)

**Power Chairs**

**Model #: GP-205**

**\$2,499**

✓ Free Delivery \*



- Full-size charcoal gray captain's seat standard
- 50 amp, PG controller standard
- Black, non-scuffing tires with silver rims
- Available in Red
- Front wheel drive with in-line motors
- Optional 20"x18" seat available

For Additional Information Please Click [Here](#)

**Alante Sport**

Weight Capacity	Seat Width/Depth	Turning Radius	Driving Range	Max Speed
300 lbs	18"/18"	25"	16.5 miles	3.5 mph

\*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00  
(mileage charges may apply-see above)

**Power Chairs**

**Model #: GP-605**

**\$3,399**

✓ Free Delivery \*



- Desk-length flip-up armrests
- Center-wheel drive
- Fold-flat seat hinge for a compact fold
- Height and angle adjustable footrests
- Driving range of 11.25 miles
- Dynamic Shark 60 amp controller comes standard
- Two U-1 batteries
- Electronic speed control as well as electro-mechanical brakes
- Made in the U.S

For Additional Information Please Click [Here](#)

**Compass Sport**

Weight Capacity	Seat Width/Depth	Ground Clearance	Driving Range	Max Speed
300 lbs	16"-22"/16"-20"	3.2"	11 miles	4.5 mph

\*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00  
(mileage charges may apply-see above)

**Power Chairs**

**Model #: GP-620**

**\$4,599**

✓ Free Delivery \*



- Standard two-tone captain's seat with our exclusive fold-flat seat hinge design
- Exclusive slide seat mechanism for easy seat removal and attachment
- Adjustable headrest and armrests
- 75 amp peak Dynamic Shark controller standard
- Off-board charger

For Additional Information Please Click [Here](#)

**Compass HD**

Weight Capacity	Seat Width/Depth	Turning Radius	Driving Range	Max Speed
450 lbs	20"/20"	20"	12 miles	4.5 mph

\*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00  
(mileage charges may apply-see above)