

Bathroom Safety

Model Number: RTL12202KDR

\$55

✓ Free Delivery *

Bathroom Safety Shower Tub Bench Chair



- Aluminum frame is lightweight, durable and corrosion proof
- Angled legs with suction style tips provide additional stability
- Drainage holes in seat reduce slipping
- Easy to clean
- Easy, tool free assembly of seat and legs
- Environment friendly product
- Pinch free cover allows for push pins to be depressed without pinching fingers; the audible "click" will be heard when secure
- Support collar prevents leg movement

For Additional Information Please Click [Here](#)

Outside Legs	Seat Dimensions	Seat Height	Weight	Weight Limit
19.5"(W)x19.5"(D)	19.25"(W)x11.5"(D)	15.5"-19.5"	7 lbs	400 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: 12445KD-2

\$95 Each

✓ Free Delivery *

Bath Bench with Padded Arms



- Aluminum frame is lightweight, durable and corrosion proof
- Blow molded bench provides comfort and strength
- Drainage holes in bench reduce slipping
- Legs are height adjustable to accommodate large and small individuals
- Removable padded armrests provide comfort and added security
- Tool free assembly

For Additional Information Please Click [Here](#)

Outside Legs	Seat Width	Seat Depth	Seat Height	Weight Limit
17.75"(W)x19.25"(D)	16"	14.5",16.25"	16"-21"	300 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: RTL12203KDR

\$39

✓ Free Delivery *

Bathroom Safety Shower Tub Bench Chair



- Aluminum frame is lightweight, durable and corrosion proof
- Angled legs with suction style tips provide additional stability
- Drainage holes in seat reduce slipping
- Easy to clean
- Easy, tool free assembly of seat and legs
- Environment friendly product
- Pinch free cover allows for push pins to be depressed without pinching fingers; the audible "click" will be heard when secure
- Support collar prevents leg movement

For Additional Information Please Click [Here](#)

Outside Legs	Seat Dimensions	Seat Height	Weight	Weight Limit
19.5"(W)x19.5"(D)	19.25"(W)x11.5"(D)	15.5"-19.5"	5 lbs	400 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: RTL12004KD

\$39

✓ Free Delivery *

Adjustable Height Bath Stool



- Aluminum frame is lightweight, durable and corrosion proof
- Impact resistant, composite seat is crack-proof and tarnish resistant
- Support collar reduces rattle

For Additional Information Please Click [Here](#)

Outside Legs	Seat Dimensions	Seat Height	Weight	Weight Limit
12.5"(W)x12.5"(D)	13"(W)x13"(D)	13.5"-21"	3 lbs	300 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: 12037

\$45

✓ Free Delivery *

Handheld Shower Head Spray with Diverter Valve



- Comes with diverter valve that allows either the regular shower head to operate or be used as a hand held shower spray
- On / Off switch built into handle for easy access
- Wall holder and extra-long 80" white reinforced nylon hose provides extra convenience

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: 12011KD-1

\$95

✓ Free Delivery *

Plastic Transfer Bench with Adjustable Backrest



- All transfer benches are reversible to accommodate any bathroom
- Durable blow-molded plastic bench and backrest
- Extra large, locking suction cups provide added safety
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs
- New "A" frame construction provides additional stability
- No exposed hardware
- Tool-Free Assembly requires connecting Back, Legs and Arm

For Additional Information Please Click [Here](#)

Outside Legs	Seat Dimensions	Seat Height	Weight	Weight Limit
32"(W)x23"(D)	26"(W)x18.5"-19.5"(D)	17.5"-22.5"	10 lbs	400 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: RTL12075

\$175

✓ Free Delivery *

Folding Universal Sliding Transfer Bench



- Comes with safety net under seat to prevent the dropping of soap, shampoo, etc
- Folds flat for transport and storage
- Seat belt for added safety
- Standard removable soap dish
- Transfer Bench slides left and right along frame for easy entry into and out of tub

For Additional Information Please Click [Here](#)

Seat Dimensions	Height	Weight	Weight Limit
22.5"(W)x18.5"(D)	21"-25"	16 lbs	300 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: 12025KD-1

\$199

✓ Free Delivery *

Heavy Duty Bariatric Plastic Seat Transfer Bench



- 1" Aluminum frame is lightweight, sturdy and corrosion resistant
- 500 lb Weight Capacity
- Backrest is reversible to accommodate any bathroom
- Durable blow-molded plastic bench and backrest
- Extra large locking suction cups provide added safety
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs
- Tool Free Assembly

For Additional Information Please Click [Here](#)

Outside Legs	Seat Dimensions	Seat Height	Weight	Weight Limit
32.5"(W)x24"(D)	26.5"(W)x16.5"(D)	18.5"-23.5"	14 lbs	500 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: BL200-DR

\$985

✓ Free Delivery *

AquaJoy Saver Bathlift



- 300 lb weight capacity
- Easy release suction cups on base provides maximum stability in the tub
- Fits almost any style bath
- Fits most style bathtubs
- Fully waterproof, floating battery hand control with easy to use convex and concave operating buttons
- Makes getting in and out of the bathtub safe, simple, and worry-free
- Non-reclining, sturdy one piece design maximizes legroom
- Safe, will not lower unless there is sufficient power to raise the user back out of the tub

For Additional Information Please Click [Here](#)

Seat Height	Seat Width with Flap	Weight	Weight Limit
2.7"-18.3"	28.5"	24.2 lbs	350 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: BL100-DR

\$1,625

✓ Free Delivery *

AquaJoy Premier Plus Reclining Bathlift



- Comfortable, machine washable fabric covers
- Easy release suction cups on base provides maximum stability in the tub
- Fits most style bathtubs
- Fully waterproof, floating battery hand control with easy to use convex and concave operating buttons
- Heavy duty seat and base providing maximum stability and support
- Makes getting in and out of the bathtub safe, simple and worry-free
- Many convenient options available
- Safe, will not lower unless there is sufficient power to raise the user back out of the tub
- Simple to understand charging process

For Additional Information Please Click [Here](#)

Seat Height Range	Seat Width with Flap	Weight	Weight Limit
3.1"-18.1"	27"	28.9 lbs	375 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: 12014

\$60

✓ Free Delivery *

Raised Toilet Seat



- Designed for individuals who have difficulty sitting down or standing up from the toilet
- Fits most standard and elongated toilets
- Lightweight and portable
- Standing locking mechanism ensures safety
- Wide opening in front and back for personal hygiene

For Additional Information Please Click [Here](#)

Seat Dimensions	Weight	Weight Limit
17"(W)x16.5"(D)x4.5"(H)	4 lbs	300 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: 12008KDR

\$82

✓ Free Delivery *

Raised Toilet Seat with Removable Padded Arms



- Arms can be removed or added as needed
- Clamping mechanism ensures secure locking onto toilet
- Heavy duty molded plastic construction provides additional strength and durability
- Padded arms can be easily removed without tools
- Prevents shifting of seat while transferring

For Additional Information Please Click [Here](#)

Seat Dimensions	Width Between Arms	Weight	Weight Limit
15"(W)x16.5"(D)x5"(H)	18.5"	5 lbs	300 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: 12001KD-1

\$58

✓ Free Delivery *

Toilet Safety Frame



- Anodized aluminum is sturdy and lightweight
- Arms are height and width adjustable
- Powder coated aluminum bracket easily attaches frame to toilet
- Waterfall armrests provide additional comfort and support

For Additional Information Please Click [Here](#)

Weight Limit
300 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Model Number: 12039

\$51

✓ Free Delivery *

Bath Shower Tub Grab Rail



- Attractive, ergonomically designed tub rail provides a comfortable, safe grip
- Tool free, easy-to-install
- Width adjusts from 2.75" - 6.5"

For Additional Information Please Click [Here](#)

Adjustable Clamp	Height	Weight	Weight Limit
2.75"(smallest)/6.25"(widest)	12"	6 lbs	300 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Bathroom Safety

Chrome Knurled Grab Bar



- Knurled grab bar diameter provides a no-slip grip
- Mounts vertically or horizontally
- 1 ½" distance from wall
- Attractive retail packaging

❖ For 12 inch Grab Bar **(Model Number: 12112-3)**

\$20 Each

❖ For 16 inch Grab Bar **(Model Number: 12116-3)**

\$30 Each

❖ For 18 inch Grab Bar **(Model Number: RTL12118)**

\$35 Each

❖ For 24 inch Grab Bar **(Model Number: 12124-3)**

\$40 Each

❖ For 32 inch Grab Bar **(Model Number: 12124-3)**

\$45 Each

Free Delivery *

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.