

Commodes

Model Number: 11148KDR

\$75

✓ Free Delivery *

Folding Bedside Commode with Bucket and Splash Guard



- Blue powder coated steel welded construction increases strength and durability
- Comes complete with 12 qt commode bucket with carry handle, cover and splash shield
- Durable plastic snap on seat and lid
- Easily opens and folds
- Elastic cord allows legs to be folded up to reduce size for storage and transport
- Folds flat for convenient storage and transportation
- Legs fold up for smallest folded commode
- Plastic arms for added comfort

For Additional Information Please Click [Here](#)

Width	Seat Dimensions	Height	Weight	Weight Limit
18" Between Arms/22.25" Outside	13"(W)x14.5"(D)	15.75"-21.5"	13 lbs	350 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Commodes

Model Number: 11112-2

\$206 Each

✓ Free Delivery *

Aluminum Commode with Wheels



- Brass push button allows for easy, tool-free assembly and disassembly of back
- Non-skid, rust-resistant, swivel casters (2 with locks)
- Aluminum construction is lightweight and durable
- Angled back for support and comfort
- Full length plastic arms for added comfort
- Comes with a 7.5 qt commode bucket, splash shield and cover. Splash shield acts as pail holder

For Additional Information Please Click [Here](#)

Width	Seat Dimensions	Height	Weight	Weight Limit
19.5" Between Arms/23.5" Outside	13.75"(W)x16.5"(D)	19.25"	12 lbs	300 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Commodes

Model Number: 11125KD-1

\$202

✓ Free Delivery *

Steel Drop Arm Bedside Commode with Padded Arms



- Comes with 12 qt commode bucket with carry handle and splash shield
- Easy to assemble frame
- Easy to clean grey powder coated steel finish
- Easy-to-release drop arm mechanism allows for safe lateral patient transfers to and from commode
- Ideal for those with limited dexterity
- Legs are height adjustable
- Plastic armrests and back provide extra comfort and support

For Additional Information Please Click [Here](#)

Width	Seat Dimensions	Height	Weight	Weight Limit
19.5" Between Arms/23.5" Outside	13.75"(W)x16.5"(D)	20"-24"	18 lbs	300 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Commodes

Model Number: 11125PSKD-1

\$245

✓ Free Delivery *

Steel Drop Arm Bedside Commode with Padded Seat & Arms



- Comes with 12 qt commode bucket with carry handle and splash shield
- Easy to assemble frame
- Easy to clean grey powder coated steel finish
- Easy-to-release drop arm mechanism allows for safe lateral patient transfers to and from commode
- Ideal for those with limited dexterity
- Legs are height adjustable
- Padded open front vinyl toilet seat
- Plastic armrests and back provide extra comfort and support

For Additional Information Please Click [Here](#)

Width	Seat Dimensions	Height	Weight	Weight Limit
19" Between Arms/23.75" Outside	15.25"(W)x16.75"(D)	19"-23"	21 lbs	300 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Commodes

Model Number: 11117N-1

\$195

✓ Free Delivery *

Heavy Duty Bariatric Folding Bedside Commode Seat



- 16.5" deep snap on seat with lid
- Easily opens and folds
- Folds flat for convenient storage and transportation
- Plastic armrests provide comfort and stability

For Additional Information Please Click [Here](#)

Width	Seat Dimensions	Height	Weight	Weight Limit
24" Between Arms/31.75" Outside	13.75"(W)x16.5"(D)	15.5"-22"	19 lbs	650 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Commodes

Model Number: 11117NPS-2

\$215 Each

✓ Free Delivery *

Heavy Duty Bariatric Folding Bedside Commode Seat



- 16.5" deep snap on seat with lid
- Easily opens and folds
- Folds flat for convenient storage and transportation
- Plastic armrests provide comfort and stability

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Commodes

Model Number: 11130-2

\$285 Each

✓ Free Delivery *

Deluxe Bariatric Commode



- Large, plastic backrest provides patient comfort
- Durable, heavy-duty, gray powder-coated steel tubing is sturdy and easy- to-maintain
- Large, durable snap-on seat
- Plastic armrests provide comfort and stability
- Pail easily slides in from the front of the commode
- Comes complete with 12 qt commode bucket and cover that removes from behind
- Weight Capacity 650 lbs

For Additional Information Please Click [Here](#)

Width	Seat Dimensions	Height	Weight	Weight Limit
23" Between Arms/27" Outside	23.25"(W)x18.5"(D)	15.5"-22"	17 lbs	650 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Commodes

Model Number: 11135-1

\$385

✓ Free Delivery *

Deluxe Bariatric Drop-Arm Commode



- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode
- Durable, heavy-duty, gray powder-coated steel tubing is sturdy and easy-to-maintain
- Large, durable snap-on seat
- Padded armrests provide comfort and stability
- Pail easily slides in from the front of the commode
- Comes complete with 12 qt commode bucket and cover that removes from front
- Weight Capacity 1000 lbs

For Additional Information Please Click [Here](#)

Width	Seat Dimensions	Height	Weight	Weight Limit
25.75" Between Arms/26.25" Outside	23.25"(W)x18.5"(D)	17.5"-22"	21 lbs	1000 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Commodes

Model Number: 11132-1

\$850

✓ Free Delivery *

Bariatric Drop Arm Bedside Commode Seat



- 1000 lbs Weight Capacity
- Comes complete with extra large commode bucket and cover that removes from behind
- Easy-to-release drop arm mechanism allows for safe lateral transfer to and from commode
- Extra large opening in toilet seat, 8.5" x 12"
- Extra large, padded deep seat provides comfort
- Heavy duty, white powder coated steel construction
- Padded armrests provide comfort and support

For Additional Information Please Click [Here](#)

Width	Seat Dimensions	Height	Weight	Weight Limit
36.6" Between Arms/40" Outside	36"(W)x21"(D)	NA-18.5"	49 lbs	1000 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.