

Hospital Beds

Model Number: 15624

\$275

✓ Free Delivery *

Trapeze



For Additional Information Please Click [Here](#)

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Hospital Beds

Ultra Light Electric Beds



- Channel frame construction provides superior strength and reduced weight
- In the event of a power failure, a 9 volt battery, located in the motor, can lower the head and foot sections nine times
- Labels and color-coded springs assure accurate installation of side rails for patient safety
- The motor assembly can be installed and removed with patient in bed
- The wood grain panels are break and scratch resistant
- Unique motor is completely self-contained to reduce weight and noise

| | |
|--|----------------|
| ❖ Full Rails & Fully Electric (Model Number: 15033BV-FR) | \$1,000 |
| ❖ Half Rails & Fully Electric (Model Number: 15033BV-HR) | \$1,000 |
| ❖ Bed, Rails, Mattress & Fully Electric (Model Number: 15033BV-PKG) | \$1,265 |
| ❖ Full Rails & Semi Electric (Model Number: 15030BV-FR) | \$950 |
| ❖ Half Rails & Semi Electric (Model Number: 15030BV-HR) | \$950 |
| ❖ Bed, Rails, Mattress & Semi Electric (Model Number: 15030BV-PKG) | \$1,215 |
| ❖ Bed, Half Rails, Mattress & Semi Electric (Model Number: 15030BV-PKG-1) | \$1,215 |

For Additional Information Please Click [Here](#)

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Hospital Beds

Model Number: 15300

\$2,250

✓ Free Delivery *

Full Electric Bariatric Hospital Bed



- All steel construction
- Emergency manual crank included
- Hand control (included) provides for multiple bed positioning for patients
- Heavy duty frame ensures strength and patient safety
- Larger sleep surface than conventional bed
- Low air loss/alternating pressure mattress available
- No tools necessary for set-up
- Quiet, smooth operation
- Split pan design with removable bed ends is easy to set-up

For Additional Information Please Click [Here](#)

| Bed Height | Unit Dimensions | Sleeping Surface | Weight Limit |
|------------|-----------------|------------------|--------------|
| 18.25" | 42"(W)x88"(L) | 80" | 600 lbs |

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Hospital Beds

Model Number: 15301

\$700

✓ Free Delivery *

Bariatric Foam Mattress 42"



- High-Density foam mattress for the bariatric patient
- Meets Federal Fire Code CFR 16 part 1633
- The nylon cover is fluid-resistant

For Additional Information Please Click [Here](#)

| |
|---------------------|
| Weight Limit |
| 585 lbs |

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Hospital Beds

Model Number: 15302

\$4,000

✓ Free Delivery *

Full Electric Heavy Duty Bariatric Hospital Bed



- All steel construction
- Emergency manual crank included
- Hand control (included) provides for multiple bed positioning for patients
- Larger sleep surface than conventional bed
- Low air loss/alternating pressure mattress available
- Split pan design with removable bed ends is easy to set-up

For Additional Information Please Click [Here](#)

| Bed Height | Unit Dimensions | Weight | Weight Limit |
|---------------|-----------------|---------|--------------|
| 18.25"-26.25" | 48"(W)x88"(L) | 312 lbs | 750 lbs |

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Hospital Beds

Model Number: 15303

\$4,500

✓ Free Delivery *

Full Electric Super Heavy Duty Bariatric Hospital Bed



- All steel construction
- Emergency manual crank included
- Hand control (included) provides for multiple bed positioning for patients
- Larger sleep surface than conventional bed
- Low air loss/alternating pressure mattress available
- Split pan design with removable bed ends is easy to set-up

For Additional Information Please Click [Here](#)

| Bed Height | Unit Dimensions | Weight | Weight Limit |
|---------------|-----------------|---------|--------------|
| 18.25"-26.25" | 54"(W)x88"(L) | 332 lbs | 1,000 lbs |

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Hospital Beds

Model Number: 15310

\$1,000

✓ Free Delivery *

Bariatric Foam Mattress 48"



- High-Density foam mattress for the bariatric patient
- Meets Federal Fire Code CFR 16 part 1633
- The nylon cover is fluid-resistant

For Additional Information Please Click [Here](#)

| Weight Limit |
|--------------|
| 665 lbs |

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Hospital Beds

Model Number: 15312

\$1,150

✓ Free Delivery *

Bariatric Foam Mattress 54"



- High-Density foam mattress for the bariatric patient
- Meets Federal Fire Code CFR 16 part 1633
- The nylon cover is fluid-resistant

For Additional Information Please Click [Here](#)

| Weight Limit |
|--------------|
| 750 lbs |

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Hospital Beds

Model Number: 15624

\$275

✓ Free Delivery *

Trapeze



For Additional Information Please Click [Here](#)

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Hospital Beds

Model Number: 15201BV

\$140

✓ Free Delivery *

Half-Length bed rail



- Unique design prevents individuals from being lodged between bars
- Ideal for maximum protection when bed is elevated
- Fits most hospital beds with steel frame
- Attractive brown-vein finish is easy-to-maintain and less institutional looking
- Constructed of 1" steel
- Not for use with "Pan-Style" bed surface
- Easy-to-install

For Additional Information Please Click [Here](#)

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Hospital Beds

Model Number: 15208BV

\$185

✓ Free Delivery *

Tool Free Adjustable Half Length Bed Rail



- Attractive brown vein finish
- Constructed of 1" steel
- Easy to install
- Fits most beds with steel frame
- Half rail was designed to comply as part of the FDA Hospital Bed System Dimensional and Assessment Guidance to Reduce Entrapment
- Ideal for protection when bed is elevated
- Not for use with "Pan-Style" bed surfaces
- Tool free width adjustment allows the rail to be adjusted in 2" increments from 36" to 42"

For Additional Information Please Click [Here](#)

| Length Adjustment |
|-------------------|
| 36",38",40",42" |

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)


Hospital Beds

Model Number: 16500BV

\$160

✓ Free Delivery *

Home Bed Style Adjustable Length Bed Rails



- Adjusts between 36" and 72"
- Attractive brown vein finish
- Constructed of 1" steel
- Crossbars are easy to install and remove
- Crossbars come with extensions that will accommodate a KING SIZE BED
- Easily adjusts up or down with spring-loaded release
- Fits all home style beds
- Rail length can be adjusted without tools to 37", 43" or 57" long

For Additional Information Please Click [Here](#)

| Rail Length Adjustment | Crossbar Adjustment |
|------------------------|---------------------|
| 37", 43", or 57" | 36"-72" |

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)