

Patient Room

Model Number: 17100-BV

\$415

✓ Free Delivery *

Hip High Chair with Padded Seat



- Attractive vinyl seat and backrest is soft, comfortable and water resistant
- Attractive, durable brown vein finish
- Contoured armrests are comfortable and durable
- Designed for post-hip surgery residents

For Additional Information Please Click [Here](#)

Width Between Arms	Overall Dimensions	Weight	Weight Limit
17.5"	22"(W)x40.75"(H)x29"(D)	37 lbs	300 lbs

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Patient Room

Model Number: 13023SV

\$795

✓ Free Delivery *

Hydraulic Patient Lift with Six Point Cradle



- 6 point cradle design
- Adjustable width base
- Easy to operate caster brakes provide additional security
- High performance hydraulics raise or lower individuals gradually and safely
- Includes sling chains
- Safely raises or lowers individuals from any stationary position
- Steel construction provides maximum strength

For Additional Information Please Click [Here](#)

Base Width	Base Height	Weight	Weight Limit
42"(open)/22.5"(closed)	6.5"	72 lbs	450 lbs

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$100.00
(mileage charges may apply-see above)

Model Number: 13240

\$2,545

✓ Free Delivery *

Electric Patient Lift with Rechargeable Battery (with 6-Point Cradle)



- Audible warning when battery power falls below 50%
- Battery does not need To charge, simply plug in unit
- Battery does not need while plugged in
- Foldable mast and base for easy storage and transportation
- For use with 4 point and 2 point slings
- Heavy gauge construction
- Ideal for Homecare use
- Jumbo actuator pump with emergency button can switch to manual mode lowers patient down safely
- The "CHARGE" light will automatically shut off when battery is fully charged

For Additional Information Please Click [Here](#)

Base Width	Base Height	Weight	Weight Limit
41"(open)/24"(closed)	4.75"	108 lbs	450 lbs

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$100.00
(mileage charges may apply-see above)

Patient Room

Model Number: 13242

\$2,720

✓ Free Delivery *

Electric Patient Lift with Rechargeable Battery (with Rechargeable, Removable Battery)



- Audible warning when battery power falls below 50%
- Battery does not need To charge, simply plug in unit
- Battery does not need while plugged in
- Foldable mast and base for easy storage and transportation
- For use with 4 point and 2 point slings
- Heavy gauge construction
- Ideal for Homecare use
- Jumbo actuator pump with emergency button can switch to manual mode lowers patient down safely
- The "CHARGE" light will automatically shut off when battery is fully charged

For Additional Information Please Click [Here](#)

Base Width	Base Height	Weight	Weight Limit
41"(open)/24"(closed)	4.5"	105 lbs	450 lbs

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$100.00

(mileage charges may apply-see above)

Patient Room

Model Number: 13244

\$3,770

✓ Free Delivery *

Bariatric Electric Patient Lift (with Rechargeable Battery and Four Point Cradle)



- Base legs adjust easily and lock securely into open position with padded shifter handle
- Dual, 3" heavy duty front casters
- Extra wide coated swivel bar with 360 degrees rotation with 4 point hook up
- High lift range allows transfer to high surfaces
- Ideal for homecare use
- Jumbo actuator pump with emergency button can switch to manual mode lowers patient down safely
- Lift will not operate while plugged in
- To charge, simply plug in unit

For Additional Information Please Click [Here](#)

Base Width	Base Height	Weight	Weight Limit
45.75" (open)/25" (closed)	4.5"	132 lbs	600 lbs

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$100.00
(mileage charges may apply-see above)

Patient Room

Model Number: 13245

\$4,895

✓ Free Delivery *

Bariatric Electric Patient Lift (with Rechargeable Battery and Four Point Cradle)



- Base legs adjust easily and lock securely into open position with padded shifter handle
- Dual, 3" heavy duty front casters
- Extra wide coated swivel bar with 360 degrees rotation with 4 point hook up
- High lift range allows transfer to high surfaces
- Ideal for homecare use
- Jumbo actuator pump with emergency button can switch to manual mode lowers patient down safely
- Lift will not operate while plugged in

For Additional Information Please Click [Here](#)

Base Width	Base Height	Weight	Weight Limit
45.75" (open)/25" (closed)	4"	129 lbs	600 lbs

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$100.00
(mileage charges may apply-see above)