

**Rollators**

**Model Number: 10257RD-1/10257BL-1**

**\$159**

✓ Free Delivery \*

**Durable 4 Wheel Rollator with 7.5" Casters**



- Comes with durable, padded seat
- Removable, hinged, padded backrest can be folded up and down
- Comes with a basket to store personal items
- Durable steel frame
- 7.5" non-marring casters are ideal for indoor and outdoor use
- Easy-to-use deluxe loop locks
- Lever allows for independent handle height and angle adjustment

For Additional Information Please Click [Here](#)

Weight Limit	Handle Height	Seat Dimensions	Weight	Unit Dimensions
300 lbs	31"-37"	12"(L)x12"(W)	18.6 lbs	23.5"(W)

\*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

**Rollators**

**Model Number: R726BK**

**\$179**

✓ Free Delivery \*

**Rollator Walker with Fold Up and Removable Back Support and Padded Seat**



- 6" black non-marring wheels with soft grip tires are ideal for indoor and outdoor use
- Basket available for additional storage sold separately (item # 9501S72606)
- Comes with new seamless padded seat with zippered pouch under seat
- Comes with standard carry pouch under the seat for added privacy and security of personal items
- Easy to use deluxe loop locks
- Ergonomic handles are easy to grip and are adjustable in height
- Removable, hinged, padded backrest can be folded up or down
- Serrated brakes
- Color options: blue, red, green, and silver.

For Additional Information Please Click [Here](#)

Weight Limit	Handle Height	Seat Dimensions	Weight	Unit Dimensions
300 lbs	32" - 37"	14" (W) x 20" (H) x 12" (D)	13 lbs.	24.5" (W) x 27.5" (L)

\*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

## Rollators

**Model Number: R728BK**

**\$196**

✓ Free Delivery \*

### Rollator Walker with Fold Up and Removable Back Support and Padded Seat



- Comes standard with pouch
- Comes with padded seat with zippered pouch under seat
- Easy to use loop locks
- Handles are adjustable in height
- Large 7.5" casters are ideal for indoor and outdoor use
- Optional basket (item # 9501S72606)
- Removable, hinged, padded backrest can be folded up or down as necessary
- Color options: blue and red

For Additional Information Please Click [Here](#)

Weight Limit	Handle Height	Seat Dimensions	Weight	Unit Dimensions
300 lbs	33"-38"	14"(W) x 23"(H) x 12"(D)	14 lbs.	24"(W) x 28"(L)

\*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

## Rollators

Model Number: 510

**\$206**

✓ Free Delivery \*

### Mimi Lite Push Brake Rollator Walker



- 6" casters with soft-grip tires are ideal for indoor and outdoor use
- Comes standard with carry pouch
- Comfortable padded seat with zippered inner pouch
- Easy, one hand folding
- Handles are adjustable in height
- Optional basket available (Item #805)
- Push down lock operates easily and ensures safety
- Tool Free removable padded back rest for comfort

For Additional Information Please Click [Here](#)

Weight Limit	Handle Height	Seat Dimensions	Weight	Unit Dimensions
300 lbs	32" - 36"	12" (W) x 22" (H) x 13" (D)	15 lbs	24" (W) x 25" (L)

\*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

**Rollators**

**Model Number: 795B**

**\$355**

✓ Free Delivery \*

**Duet Transport Wheelchair Rollator Walker**



- A safety lock release is used to store the arms and a hook & loop fastener strap is used to store the footrest while using the Duet for ambulation
- An individual can ambulate independently or be pushed safely by a care giver
- Combines the features of a rollator or transport chair in one unit
- Comfortable padded seat with a large under seat pouch for safe storage of personal items
- Large 8" casters are used for both indoor and outdoor activities
- The Duet has locking flip up armrests and a flip down footrest to aid the user while the Duet is used in the transport chair mode
- The folding frame has a lock out feature to insure the unit stays open while in use and insures the unit will stay folded during transport

For Additional Information Please Click [Here](#)

Weight Limit	Handle Height	Seat Dimensions	Weight	Unit Dimensions
300 lbs	31.5" - 37"	13.5" (W) x 13" (D)	20 lbs	25.5" (W) x 37" (H) x 26" (L)

\*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

**Rollators**

**Model Number: 10289BL**

**\$186**

✓ Free Delivery \*

**3 Wheel Rollator Walker with Basket Tray and Pouch**



- Adjustable handle height
- Comes standard with basket, tray and pouch
- Comes with frame folded and compact
- Easy one hand folding
- Lightweight aluminum frame
- Lightweight, solid 8" wheels for indoor or outdoor use
- Special loop lock made of internal aluminum casting operates easily and ensures safety

For Additional Information Please Click [Here](#)

Weight Limit	Handle Height	Weight	Unit Dimensions
300 lbs	31" - 38"	11 lbs.	25" (W) x 24" (L)

\*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

## Rollators

Model Number: 10252BL

**\$295**

✓ Free Delivery \*

### Bariatric Rollator with 7.5" Wheels



- Black 7 5" casters with soft-grip tires are ideal for indoor and outdoor use
- Soft padded oversized seat
- Comes with standard carry pouch under the seat for added privacy and security of personal items
- Tool-free removable padded backrest for comfort
- Special loop lock made of internal aluminum casting operates easily and ensures safety **(Figure A)**
- Handles easily adjustable in height with self-threading knob **(Figure B)**
- Limited Lifetime Warranty

For Additional Information Please Click [Here](#)

Weight Limit	Handle Height	Seat Dimensions	Weight	Unit Dimensions
500 lbs	35" - 39"	18" (W) x 14" (D)	26 lbs.	30" (W) x 24.5" (L)

\*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

**Rollators**

**Model Number: 10216RD-1**

**\$290**

✓ Free Delivery \*

**Bariatric Rollator with 8" Wheels**



- 8" casters with soft-grip tires are ideal for indoor and outdoor use
- Soft padded oversized seat
- Comes with pouch under seat for added privacy and security of personal items
- Tool-free removable padded backrest for comfort
- Special loop lock made of internal aluminum casting operates easily and ensures safety
- Handles easily adjustable in height with self-threading knob
- Limited Lifetime Warranty

For Additional Information Please Click [Here](#)

Weight Limit	Handle Height	Seat Dimensions	Weight	Unit Dimensions
400 lbs	34" - 39"	19" (W) x 13" (D)	23.5 lbs.	28" (W) x 27.5" (L)

\*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.