

Electrotherapy

Model Number: AGF-3E

\$68

✓ Free Delivery *

Portable Dual Channel TENS Unit with Electrodes and Carry Case



- Comes complete with carrying case, 2 lead wires, 9-volt battery, 1 package of 4 pre-gelled electrodes and operating manual
- Designed for rehabilitation and physical therapy
- Dual isolated channels
- For chronic, acute and postoperative pain
- Intensity control is adjustable from 0-80 mA, 500 load
- Metal belt clip
- Pulse frequency adjustable from 2-150 Hz
- Pulse width is adjustable from 40-250 micro-seconds
- Wave form: asymmetrical, bi-phase, square pulse

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Electrotherapy

Model Number: AGF-3X

\$85

✓ Free Delivery *

Portable Dual Channel TENS Unit with Electrodes and Carry Case



- Adjustable 0-80 mA, 500 ohm load
- Comes complete with carrying case, 2 lead wires, with deluxe connector, 9-volt battery, 1 package of 4 pre-gelled electrodes and operating manual
- Dual isolated channels
- Modulation selector offers wave delivery in normal, burst or modulation mode
- New and improved lead wire connectors will not disengage from TENS unit
- Pulse frequency adjustable from 2-150 Hz
- Pulse width is adjustable from 60-250 micro-seconds
- Wave form: asymmetrical, bi-phase, square pulse

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Electrotherapy

Model Number: AGF-5X

\$125

✓ Free Delivery *

Portable Digital Dual Channel 5 Mode TENS Unit with Carrying Case and Electrodes



- Comes complete with carrying case, 2 lead wires, 9-volt battery, 1 package of 4 pre-gelled electrodes and operating manual
- Digital wave technology provides a highly efficient delivery of the wave signal
- Dual isolated channels
- Intensity control is adjustable to 80 mA, 500 ohm load
- Modulation selector offers 5 modes of wave delivery, Normal, Burst and 3 Modulation modes
- Pulse frequency adjustable from 0
- Pulse width is adjustable from 30-260 micro-seconds
- Two strength duration modes & 5-150 Hz
- Wave form: Asymmetrical, bi-phase, square pulse

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Electrotherapy

Model Number: AGF-602

\$90

✓ Free Delivery *

Portable Dual Channel TENS Unit with Timer and Electrodes



- Comes complete with carrying case, 2 lead wires, 9-volt battery, 1 package of 4 pre-gelled electrodes and operating manual
- Dual isolated channels
- Modulation selector offers wave delivery in Normal, Burst or Modulation modes
- Plastic belt clip
- Pulse amplitude adjustable from 0-80 mA, 500 ohm load
- Pulse frequency adjustable from 2-150 Hz
- Pulse width is adjustable from 30-260 micro-seconds
- Wave form: asymmetrical, bi-phase, square pulse

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Electrotherapy

Model Number: AGF-6X

\$130

✓ Free Delivery *

Deluxe Digital Electronic Muscle Stimulator



- Three Stimulation Modes: Constant, Synchronous, and Alternate
- Timer: Adjustable 15, 30, 45, 60 minutes, or continuous
- Digital wave technology provides a highly efficient delivery of the wave signal
- Wave Form: Asymmetrical rectangular biphasic pulse
- Pulse frequency adjustable from 2-150 Hz
- Pulse Amplitude: Adjustable 0-100 mA into 500 ohm load each channel
- Pulse width is adjustable from 50 uS ~ 300 uS, 10 uS/step
- Comes complete with carrying case, 2 lead wires, 9-volt battery, 1 package of 4 pre-gelled electrodes and operating manual

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Electrotherapy

Model Number: AMS-4N

\$120

✓ Free Delivery *

Portable EMS with Timer and Carrying Case



- Comes complete with carrying case, 2 lead wires, 9-volt battery, 1 package of 4 pre-gelled electrodes and operating manual
- Contraction time is adjustable from 2-30 seconds
- Dual isolated channels
- Intensity control is adjustable from 0-80 mA, 500 ohm load
- Metal belt clip
- Pulse rate is adjustable from 2-150 Hz
- Pulse width is selectable: 100, 200, or 300 micro-seconds
- Ramp time is adjustable from 1-6 seconds
- Relaxation time is adjustable from 1-30 seconds

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.