

Respiratory

Model Number: 18070

\$45

✓ Free Delivery *

Pacifica Elite Nebulizer



- Retail-packaged compressor neb is a workhorse with higher liter flows, shorter treatment times and quieter operation than similar economy compressor nebs
- Ultra compact design makes it extremely portable
- Powerful piston pump is efficient, dependable and quiet
- Delivers consistent particle size of <5 microns and fast treatment times
- Retail packaged

For Additional Information Please Click [Here](#)

Method of Nebulization	Unit Dimensions	Weight	Max Pressure	Operating Pressure	Liter Flow
Piston Compressor	7.5"(W)x4"(H)x6.8"(D)	4.1 lbs	35 psi	= 10 psi	= 8 lpm

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Respiratory

Model Number: HUM 001

\$5

✓ Free Delivery *

Humidifier Bottle



- Drive's multi-threaded lid ensures a secure fit every-time, preventing leaks and service calls
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem: it eliminates shredded aluminum diffuser so often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents
- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle has a 500 ml water capacity
- There are finger grips at the top of the bottle, making the lid easier to screw on/off for the patient

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Respiratory

Model Number: MASK 002A [Adult Oxygen Mask]

\$5

✓ Free Delivery *



Model Number: MASK 002P [Pediatric Oxygen Mask]

\$5

✓ Free Delivery *



Model Number: MASK 003A-1 [Adult Oxygen Mask with 7" Tubing]

\$85

✓ Free Delivery *



Respiratory

Model Number: 525DS

\$750

✓ Free Delivery *

DeVilbiss 5 Liter Oxygen Concentrator



- Built upon a legacy of reliable, field-proven concentrator models, the DeVilbiss 5 Liter Oxygen Concentrator was designed with enhanced quality, durability and simplicity in mind
- Patented DeVilbiss Turn-Down Technology minimizes wear on internal components, reduces power consumption by 15% or more and extends the life of the concentrator
- Exclusive DeVilbiss OSD (Oxygen Sensing Device), standard on every unit, ensures patient safety and reliability for longer service intervals
- Readily accessible patient controls, protected cannula fitting and recessed humidifier nook to prevent damage
- Alarms Audible and visual high/low pressure, low flow, low oxygen, power fail, Oxygen Sensing Device
- Sold separately, the Transfill Caddy is designed with a steel frame to simply connect the 18 x 14 inch shelf, designed to support up to 40 lbs, to the DeVilbiss 5L Concentrator.

For Additional Information Please Click [Here](#)

Power Consumption	Product Dimensions	Weight	Oxygen Percentage	Outlet Pressure	Delivery Rate
310 Watt Average, 275 Watts @ 1.2 LPM and Below	24.5"(H)x13.5"(W)x 12"(D)	36 lbs	1-5 LPM, 93% +/- 3%	8.5 psig	0.5 to 5 LPM

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Respiratory

Model Number: 525DS-Q

\$1,100

✓ Free Delivery *

[Ultra Quiet] DeVilbiss 5 Liter Oxygen Concentrator



- Built upon a legacy of reliable, field-proven concentrator models, the DeVilbiss 5 Liter Oxygen Concentrator was designed with enhanced quality, durability and simplicity in mind
- Patented DeVilbiss Turn-Down Technology minimizes wear on internal components, reduces power consumption by 15% or more and extends the life of the concentrator
- Exclusive DeVilbiss OSD (Oxygen Sensing Device), standard on every unit, ensures patient safety and reliability for longer service intervals
- Readily accessible patient controls, protected cannula fitting and recessed humidifier nook to prevent damage
- Alarms Audible and visual high/low pressure, low flow, low oxygen, power fail, Oxygen Sensing Device
- Sold separately, the Transfill Caddy is designed with a steel frame to simply connect the 18 x 14 inch shelf, designed to support up to 40 lbs, to the DeVilbiss 5L Concentrator.

For Additional Information Please Click [Here](#)

Power Consumption	Product Dimensions	Weight	Oxygen Percentage	Outlet Pressure	Delivery Rate
310 Watt Average, 275 Watts @ 1.2 LPM and Below	24.5"(H)x13.5"(W)x 12"(D)	36 lbs	1-5 LPM, 93% +/- 3%	8.5 psig	0.5 to 5 LPM

*Free curbside delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

White glove delivery, set up, and training for ONLY \$75.00
(mileage charges may apply-see above)

Respiratory

Model Number: 18707

\$95

✓ Free Delivery *

Pediatric Fingertip Pulse Oximeter



- Four-way OLED with six display modes
- Prescription required
- Portable, non-invasive spot check of oxygen saturation and pulse rate
- Uses two AAA batteries
- Low battery indicator
- Easy operation auto-off feature
- Lanyard, silicon case, and hang bag included
- Finger size range: 0.3" - 1.0"

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Respiratory

Model Number: MQ3000

\$90

✓ Free Delivery *

Fingertip Pulse Oximeter



- Large easy-to-read multi-view LED display screen
- Displays SpO2 and pulse rate
- Low power consumption, Auto-off function
- Approximately 30 hours of normal operation
- Accommodates a wide range of finger sizes from pediatric to adult
- Includes lanyard and two AAA batteries.
- Perfect for sports enthusiasts to measure oxygen saturation and pulse rate

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Respiratory

Model Number: MQ9000

\$250

✓ Free Delivery *

Gentle Feed Single Channel Breast Pump



- LED Display provides status of programs and functions
- Adjustable Multi-level Suction
- Includes Soft Silicone Cushion Breast Cup for Optimum Comfort
- Anti-backflow System Prevents Milk from Entering Motor
- Bottle Feeding Adapter Kit & Storage Container with Lid
- Bisphenol A (BPA) Free
- Convenient Carrying Bag
- May be powered with four AA Alkaline Batteries
- Includes LED light display

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Respiratory

Model Number: MQ9130

\$219

✓ Free Delivery *

Gentle Feed + Dual Channel Breast Pump



- Choice of Stimulation or Expression modes.
- Control suction and speed to suit your comfort level with a total of 28 pumping state options.
- Piston Pump System imitates baby's sucking rhythm.
- Comfortable & painless while pumping
- Silent for discreet pumping.
- Bisphenol A (BPA) Free.
- Anti-backflow system prevents milk from entering motor.
- Includes bottle feeding adapter kit & storage containers with lids.
- May be powered with four AA Alkaline batteries or the included A/C power supply.

For Additional Information Please Click [Here](#)

Power	Type	Suction Levels	Suction Speed
9V AC/4 AAs	Double	8	5 Preset Speeds

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

