

Model Number: CE 1000 NBL

\$300

✓ Free Delivery *

Adult and Bariatric Anterior Safety Rollers



- Height adjustable in 1" increments
- Width adjustable
- Adjustable-brake spring tension
- Patient is situated in the center of the Safety Roller, which provides maximum support and improved posture
- Skid-resistant, non-marring, easy-rolling neoprene wheels Ideal for indoor and outdoor use
- Folds easily and stands on its own in the folded position

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Model Number: CE 1000 B & CE 1000 BK

\$300

✓ Free Delivery *

Anterior Safety Roller



- Adjustable brake spring tension
- Applying pressure to the handlebar activates breaking mechanism
- Height adjustable in 1" increments
- Ideal for indoor and outdoor use
- Ideal for patients with limited hand function and/or limited cognizance
- Patient is situated in the center of the Safety Roller, which provides maximum support and improved posture
- Skid resistant, non-marring, easy rolling neoprene wheels
- The Safety Roller halts in step with the patient and the breaking mechanism automatically engages when patient stumbles

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Wenzelite

Model Number: CE 1000 XL

\$400

✓ Free Delivery *

Adult and Bariatric Anterior Safety Rollers



- Height adjustable in 1" increments
- Width adjustable
- Adjustable-brake spring tension
- Patient is situated in the center of the Safety Roller, which provides maximum support and improved posture
- Skid-resistant, non-marring, easy-rolling neoprene wheels Ideal for indoor and outdoor use
- Folds easily and stands on its own in the folded position

For Additional Information Please Click [Here](#)

Weight Capacity
500 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Wenzelite



Model Number: WB 1200

\$550

✓ Free Delivery *

Wallaby Pediatric Wheelchair (12")



- Standard Transit Tie Down Brackets
- WC19 Compliant when used with WB8021 for users between 50 and 80 lbs
- Height adjustable push handles
- Flip back padded desk arms
- Pelvic belt standard
- Push to lock wheel locks
- 22" Rear wheels with 6" front casters
- Swingaway footrests with tool free height adjustment and calf strap

For Additional Information Please Click [Here](#)

Seat to Floor Height	Backrest Height	Weight Capacity	Weight Capacity When in Transport
19"	14"	150 lbs	50-80 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Model Number: WB 1400

\$550

✓ Free Delivery *

Wallaby Pediatric Wheelchair (14")



- Standard Transit Tie Down Brackets
- WC19 Compliant when used with WB8021 for users between 50 and 80 lbs
- Height adjustable push handles
- Flip back padded desk arms
- Pelvic belt standard
- Push to lock wheel locks
- 22" Rear wheels with 6" front casters
- Swingaway footrests with tool free height adjustment and calf strap

For Additional Information Please Click [Here](#)

Seat to Floor Height	Backrest Height	Weight Capacity	Weight Capacity When in Transport
19"	14"	150 lbs	50-80 lbs

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.

Marx Medical
Equipment & Supplies

Wenzelite



Model Number: WB 8021

\$95

✓ Free Delivery *

Wallaby Headrest Extension



- Provides additional head support for taller users
- Width adjustable headrest extension
- Will prevent head whiplash during rear end collision on a bus

For Additional Information Please Click [Here](#)

*Free delivery within 10 miles from our location. Over this mileage there is a \$35 fuel charge every 20 miles.